

This work is licensed under a [Creative Commons Attribution-NonCommercial-ShareAlike License](https://creativecommons.org/licenses/by-nc-sa/4.0/). Your use of this material constitutes acceptance of that license and the conditions of use of materials on this site.



Copyright 2006, The Johns Hopkins University, Patrick Breyse, and Peter S. J. Lees. All rights reserved. Use of these materials permitted only in accordance with license rights granted. Materials provided "AS IS"; no representations or warranties provided. User assumes all responsibility for use, and all liability related thereto, and must independently review all materials for accuracy and efficacy. May contain materials owned by others. User is responsible for obtaining permissions for use from third parties as needed.



JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

# Air Sampling for Particulate Matter

**Patrick N. Breysse, PhD, CIH**

**Peter S.J. Lees, PhD, CIH**

**Johns Hopkins University**




JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

# Section A

## *Introduction*

# Exposure Assessment Methods

- ◆ General air-sampling methodology



# General Sampling Methodology

- ◆ Personal sampling
  - Battery powered pump worn by worker
- ◆ Sample from “breathing zone”
- ◆ Area sampling
  - Sample from fixed location

# Concentration

$$\text{Concentration} = \frac{\text{Mass}_{\text{contaminant}}}{\text{Volume}_{\text{air}}}$$

- ◆ Contaminant mass from laboratory analysis
- ◆ Air volume from product of calibrated air flow rate and sampling time



JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

## **Section B**

***Sample Pumps and  
Volume Determination***

# Personal Sampling Pumps





# Sample Volume

$$V = Q T$$


- ◆ Where  $V$  is volume, L
- ◆  $Q$  is flow rate, L/min
- ◆  $T$  is sample time, min

$$(2.5 \text{ L/min})(480 \text{ min}) = 1,200 \text{ L}$$

$$(1,200 \text{ L})(1 \text{ m}^3 / 1000 \text{ L}) = 1.2 \text{ m}^3$$



# Pump Calibration



# Frictionless Piston Flow Meter



Calibrator

Tubing

Sample collector

Pump

# Soap Bubble Flow Meter

- ◆ Measure time required for bubble to move through specified volume in burette

$$\left( \frac{1.0\text{L}}{30 \text{ sec}} \right) \left( \frac{60 \text{ sec}}{\text{min}} \right) = 2.0\text{L} / \text{min}$$




JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

# Section C

## *Sample Collection*

# Sampling for Aerosols

- ◆ Filtration is method of choice
- ◆ Filter defined by composition, diameter, and pore size
- ◆ General classes of filters used:
  - Membrane filter
  - Nuclepore filter



Source: Peter Lees

# Filter Selection

- ◆ Membrane filter types
  - Mixed cellulose ester (MCE)
  - Teflon
  - Polyvinyl chloride
- ◆ Filter type selected primarily for compatibility with analytical method



# Filter Selection


- ◆ MCE filters
  - Clears for microscopy
  - Particle counting and sizing
  - Fiber counting (asbestos)
- ◆ PVC filters
  - Do not absorb water vapor
  - Stable for gravimetric analysis

# Filter Selection

- ◆ Teflon filters
  - Low tare weight
  - Low mass gravimetric analysis

# 37 mm Filter Cassette

- ◆ Plastic cassette holds filter paper
- ◆ Filter collects “total” dust
- ◆ Filter diameter is 37 mm



Source: Peter Lees

# 37 mm Filter Cassette




# Filter Analysis

- ◆ Gravimetric (pre-weight minus post-weight)
- ◆ Specific contaminant analysis
  - Metals
  - Silica
  - Fibers
  - Polycyclic aromatic hydrocarbons

# Impingers

- ◆ Old dust methods used microscopic counting of particles (mppcf)
- ◆ Current bioaerosol method (fungal spores and bacteria)



Source: U. S. Government




JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

## **Section D**

### ***Size-Selective Samplers***

# Exposure Assessment Methods

- ◆ Size-selective air-sampling methodology




- ◆ Size selector permits collection of more biologically-relevant aerosol fraction
  - Respirable fraction, thoracic fraction, inhalable fraction, PM10, PM 2.5



# Size Selection: Impactors

- ◆ Multistage or single stage
- ◆ Multiple fractions of varying size ranges possible
  - Can be used to determine complete size distribution
- ◆ Flow rate dependent
- ◆ Commonly used for bioaerosol assessment
- ◆ Mostly area-type measurement because high flow rate required

# S-Stage Impactor Schematic




# Size Selection: Cyclones


- ◆ Cyclones are commonly used size-selective samplers
- ◆ Used to define respirable particle sampling criteria
  - Cyclone removal mimics respiratory system when operated
  - Flow rate is critical

# Cyclone Samplers

To sample  
collection filter



Inertial  
Impaction




# 10 mm Dorr-Oliver Cyclone

Filter cassette

Cyclone inlet


Cyclone body



# Personal P10 and PM2.5 Impactors

- ◆ EPA regulates particles  $<10$   $\mu\text{m}$  (PM10) and  $<2.5$   $\mu\text{m}$  (PM2.5)
- ◆ Sampled using impactors (personal and area samplers)
- ◆ Used in non-occupational sampling associated with ambient air pollution
  - Childhood asthma studies
  - Effects of PM on mortality
  - COPD Studies

# Cyclone/Impactor Performance





JOHNS HOPKINS  
BLOOMBERG  
SCHOOL of PUBLIC HEALTH

# Section E

*Direct-Reading Instruments*



# Direct-Reading Instruments

- ◆ Based on light scattering
  - Intensity of scatter light is proportional to concentration
  - Require calibration
  - Effected by humidity
- ◆ Incorporate data logger so can get concentration over time